



**DR. WEIGERT**  
Systematic Hygiene

Medical



neodisher®  
**Septo  
Active**

Manual Disinfecting Cleaning and High-level Disinfection  
of Thermally Stable and Thermally Instable Instruments.

# neodisher® Septo Active – the Granulate for Safe Manual Reprocessing of Medical Devices.

## neodisher® Septo Active

The granulate combines high-level disinfection with active cleaning.

- For the manual disinfection of thermally stable and thermally instable instruments
- For the manual cleaning with disinfecting action in immersion baths and ultrasonic baths



8 kg Bucket  
4098 76



4 x 2 kg Bucket  
4098 81



25 x 100 g Sachet  
4098 79



50 x neodisher® test strips  
9813 20

### Easy to use:

- Prepare desired water volume (max. lukewarm)
- Add granulate according to dosing table, stir until completely dissolved
- After 15 min the working solution is ready to be used
- To determine the concentration **neodisher® test strips** are available
- 1.0 % working solution (10 g/l) contains 0.15 % peracetic acid

### Application recommendation (20 °C):

#### Disinfecting Cleaning

bactericidal, yeasticidal, active against enveloped viruses (incl. e.g. HBV, HIV, HCV)

10 g/l (1.0 %),  
5 min

#### Disinfecting Cleaning with Sporicidal Activity

bactericidal, yeasticidal, active against enveloped viruses (incl. e.g. HBV, HIV, HCV), sporicidal against *C. difficile*

10 g/l (1.0 %),  
15 min

#### Disinfection

bactericidal, mycobactericidal, fungicidal, virucidal (incl. e.g. rota viruses, noro viruses, HAV), sporicidal against *C. difficile*

20 g/l (2.0 %),  
15 min

incl. sporicidal

20 g/l (2.0 %),  
15 min\*

\* When routinely reprocessed with the immersion method with a contact time of 60 minutes material alterations, especially with plastics, cannot be excluded and can therefore shorten the service life of medical devices.

## Our Competence – Your Benefit:

- extensive microbiological activity
- optimum cleaning performance
- free of aldehydes, amines and quaternary ammonium compounds
- coated granulate – available in convenient 100 g sachets
- excellent material compatibility



- safe disinfection, tested and proved according to European standards and DGHM<sup>1</sup> methods
- no protein fixation
- based on peracetic acid, no incompatibilities with other active substances
- low dust, easily soluble – clean, precise and simple dosing
- suitable for materials such as stainless steel, anodised aluminium, plastic (incl. silicone)

<sup>1</sup>Deutsche Gesellschaft für Hygiene und Mikrobiologie (German Society of Hygiene and Microbiology)