**DR. WEIGERT**  
Systematic Hygiene

Medical

neodisher®  
**Septo Fin**

Disinfectant for the manual reprocessing of  
thermostable und thermolabile instruments.

# neodisher® Septo Fin – safe disinfection in manual instrument reprocessing.

## neodisher® Septo Fin

Formaldehyde and glyoxal free disinfectant for the manual final disinfection, excellent material compatibility, usable at all degrees of water hardness

- For the disinfection of thermostable and thermolabile instruments, including flexible endoscopes, TEE probes, MIS instruments and anaesthesia equipment in immersion baths and semi-automatic reprocessors

5 L canister  
4207 35

6 x 2 L bottle  
4207 42

### Application and dosage:

- **neodisher® Septo Fin** can be used in immersion baths and semi-automatic reprocessors. Depending on the desired activity a solution is prepared according to the application recommendations. Immerse instruments resp. endoscopes in solution and adhere to the respective contact time.
- **neodisher® Septo Fin** has a bactericidal (incl. MRSA, H. pylori), mycobactericidal (incl. tuberculosis pathogen), fungicidal, virucidal and sporicidal (C. difficile) activity.
- Disinfecting active substances in 100 g: 15 g glutaral.

### Application recommendation (20 °C):

Disinfection (clean conditions)	
bactericidal, yeasticidal (VAH, EN 13727, EN 13624, EN 14561, EN 14562)	30 ml/l (3.0 %), 15 min 20 ml/l (2.0 %), 30 min 10 ml/l (1.0 %), 60 min
mycobactericidal (EN 14348, EN 14563)	50 ml/l (5.0 %), 30 min 40 ml/l (4.0 %), 60 min
fungicidal (EN 13624, EN 14562)	50 ml/l (5.0 %), 60 min
sporicidal against C. difficile (EN 17126)	40 ml/l (4.0 %), 24 hrs.
active against enveloped viruses (EN 14476, EN 17111, RKI/DVV <sup>3</sup> )	10 ml/l (1.0 %), 15 min
virucidal (EN 14476, EN 17111)	40 ml/l (4.0 %), 15 min 30 ml/l (3.0 %), 30 min

## Our Competence – Your Benefit:

- low-foam formulation
- virucidal activity in 15 min (30 ml/l, 20 °C)
- disinfecting activity tested and proved according to European standards, VAH<sup>2</sup> and (RKI/DVV)<sup>3</sup> methods
- free of formaldehyde and glyoxal, mildly perfumed
- excellent material compatibility

- suitable for use in circulation processes, e.g. semi-automatic reprocessors
- rapid microbicidal activity
- VAH<sup>2</sup>-listed, ÖGHMP<sup>4</sup>-listed, on the list of the IHO<sup>5</sup> for disinfectants
- free of CMR<sup>1</sup>-containing disinfecting active substances, very user-friendly, pleasant odour
- suitable for materials such as stainless steel, chrome nickel steel, anodised aluminium, plastics incl. silicone, rubber, porcelain, glass, non-ferrous metals

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