



DR. WEIGERT
Systematic Hygiene

Medical



neodisher®
**Septo
Active**

Manual Disinfecting Cleaning and High-level Disinfection
of Thermally Stable and Thermally Instable Instruments.

neodisher® Septo Active – the Granulate for Safe Manual Reprocessing of Medical Devices.

neodisher® Septo Active

The granulate combines high-level disinfection with active cleaning.

- For the manual disinfection of thermally stable and thermally instable instruments
- For the manual cleaning with disinfecting action in immersion baths and ultrasonic baths



8 kg Bucket
4098 76



4 x 2 kg Bucket
4098 81



25 x 100 g Sachet
4098 79

Easy to use:

- Prepare desired water volume (max. lukewarm)
- Add granulate according to dosing table, stir until completely dissolved
- After 15 min the working solution is ready to be used
- To determine the concentration **neodisher® test strips** are available
- 1.0 % working solution (10 g/l) contains 0.15 % peracetic acid

Application recommendation (20 °C):

Disinfecting Cleaning

bactericidal, yeasticidal, active against enveloped viruses (incl. e.g. HBV, HIV, HCV)

10 g/l (1.0 %),
5 min

Disinfecting Cleaning with sporicidal activity

bactericidal, yeasticidal, active against enveloped viruses (incl. e.g. HBV, HIV, HCV), sporicidal (incl. Clostridium difficile)

10 g/l (1.0 %),
15 min

Disinfection

bactericidal, mycobactericidal, fungicidal, virucidal (incl. e.g. rota viruses, noro viruses, HAV), sporicidal (incl. Clostridium difficile)

20 g/l (2.0 %),
15 min

50 strips

50 x neodisher® test strips
9813 20

Our Competence – Your Benefit:

- extensive microbiological activity
- optimum cleaning performance
- free of aldehydes, amines and quaternary ammonium compounds
- coated granulate – available in convenient 100 g sachets
- excellent material compatibility



- safe disinfection, tested and proved according to European standards and DGHM¹ methods
- no protein fixation
- based on peracetic acid, no incompatibilities with other active substances
- low dust, easily soluble – clean, precise and simple dosing
- suitable for materials such as stainless steel, anodized aluminium, plastic (incl. silicone)

¹Deutsche Gesellschaft für Hygiene und Mikrobiologie (German Society of Hygiene and Microbiology)