



DR. WEIGERT
Systematic Hygiene

Endoscopy



neodisher
endo[®]
MED

Detergent and Disinfectant for the
Manual Pre-Treatment of Flexible Endoscopes.

neodisher endo[®] MED – manual disinfecting pre-cleaning compatible with all active substances.

neodisher endo[®] MED

Aldehyde-free and alkylamine-free detergent and disinfectant for the manual pre-treatment of flexible endoscopes; with excellent material compatibility. Compatible with disinfectants containing aldehydes and peracetic acid for subsequent reprocessing.

- For the disinfecting pre-cleaning of flexible endoscopes and endoscopic accessories in immersion baths or ultrasonic baths



5 L canister
4210 35

10 x



10 x 1 L bottle
4210 45



Excellent material compatibility – suitable for endoscopes of all leading manufacturers

Application and Dosage:

- **neodisher endo[®] MED** is used in immersion baths or ultrasonic baths. A solution is prepared in accordance with the application recommendations. For disinfecting pre-cleaning place flexible endoscopes into the solution and pre-clean according to the manufacturer's recommendations. When used in ultrasonic baths the treatment time should not exceed the time given by the instrument manufacturer. Also usable for wipe cleaning of flexible endoscopes directly after the examination (bedside cleaning) in combination with the fleece wipe dispensing system neoform[®] wipes RTF.

Active substances in 100 g:

- 18 g N,N-didecyl-N-methyl-poly(oxyethyl)ammoniumpropionate

Application Recommendation (20 °C):

Disinfecting pre-cleaning

bactericidal, yeasticidal, active against enveloped viruses (incl. e.g. HBV, HIV, CV)

10 ml/l (1.0 %),
15 min

bactericidal

(according to VAH¹, EN 13727, EN 14561), dirty conditions

10 ml/l (1.0 %),
15 min

yeasticidal

(according to VAH¹, EN 13624, EN 14562), dirty conditions

10 ml/l (1.0 %),
15 min

active against enveloped viruses

(incl. HBV, HCV, HIV) (according to RKI/DVV²), dirty conditions

10 ml/l (1.0 %),
15 min

neodisher endo[®] MED
is compatible with all common
disinfecting active substances

Our Competence – Your Benefit:

- excellent cleaning performance; does not fixate proteins
- disinfecting action confirmed in accordance with EN standards, VAH¹ and RKI/DVV² methods
- free of aldehydes and alkylamines
- without CMR³ disinfecting active substances
- excellent material compatibility



- facilitates manual cleaning; optimum pre-treatment
- increased personnel protection; improves hygiene safety when reprocessing; VAH¹ listed; included in the IHO⁴ list of virucidal products
- compatible with disinfectants containing aldehydes and peracetic acid for subsequent reprocessing steps; no layers/discolouration on endoscopes
- especially user-friendly
- suitable for stainless steel, anodized aluminium, non-ferrous metals (copper and brass), synthetic materials (incl. silicone)

¹ Verband für Angewandte Hygiene [Association of Applied Hygiene] the Control of Virus Diseases]

³ Carcinogenic, Mutagenic, toxic to Reproduction

² Testing guidelines of the Robert Koch Institute and the Deutsche Vereinigung zur Bekämpfung von Viruskrankheiten [German Association for the Control of Virus Diseases]

⁴ Industrieverband Hygiene und Oberflächenschutz [German Association for Hygiene and Surface Protection]