



**DR. WEIGERT**  
Systematic Hygiene

Medical



neodisher<sup>®</sup>  
**Septo Plus**

Detergent and Disinfectant for the Manual Reprocessing  
of Thermostable and Thermolabile Instruments.

# neodisher® Septo Plus – Efficient Cleaning and Disinfection Compatible with All Active Substances.

## neodisher® Septo Plus

Detergent and disinfectant with excellent material compatibility and compatibility with disinfectants containing aldehydes and peracetic acid for subsequent reprocessing.

- For the manual cleaning and disinfection of thermostable and thermolabile instruments
- For the disinfecting pre-cleaning in immersion baths or ultrasonic baths with optimum personnel safety



5 L Canister  
4209 35



12 x 1 L Bottle  
4209 44

### Application and Dosage:

- **neodisher® Septo Plus** is used in immersion baths or ultrasonic baths. Depending on the desired activity a solution is prepared in accordance with the application recommendations. Place instruments or endoscopes into the solution whilst observing the below-mentioned contact time
- **neodisher® Septo Plus** is bactericidal, yeasticidal and active against enveloped viruses (incl. HIV, HBV, HCV), as well as active against rotavirus and polyomavirus when used with the below-mentioned parameters
- Disinfecting active substance in 100 g: 18 g N,N-didecyl-N-methyl-poly(oxyethyl)ammoniumpropionate

### Application Recommendation (20 °C):

Disinfecting activity (dirty conditions)

|  |  |
|--|--|
| bactericidal, yeasticidal (VAH <sup>1</sup> )                                  | 10 ml/l (1.0 %), 15 min                            |
| bactericidal (EN 13727, EN 14561), yeasticidal (EN 13624, EN 14562)            | 10 ml/l (1.0 %), 5 min<br>5 ml/l (0.5 %), 15 min   |
| active against enveloped viruses (incl. HBV, HCV, HIV) (RKI/DVV <sup>2</sup> ) | 15 ml/l (1.5 %), 5 min<br>10 ml/l (1.0 %), 15 min  |
| rotavirus (EN 14476)   | 10 ml/l (1.0 %), 5 min                             |
| polyoma SV 40 (EN 14476)   | 20 ml/l (2.0 %), 15 min<br>15 ml/l (1.5 %), 30 min |

Disinfecting cleaning

|  |                         |
|--|-------------------------|
| bactericidal, yeasticidal, active against enveloped viruses (incl. e.g. HBV, HIV, HCV) | 10 ml/l (1.0 %), 15 min |
|--|-------------------------|



## Our Competence – Your Benefit:

- Disinfecting action confirmed in accordance with EN standards, VAH<sup>1</sup> and RKI/DVV<sup>2</sup> methods
- Free of aldehydes and amines
- Without CMR<sup>3</sup> disinfecting active substances, perfume-free
- Excellent cleaning performance
- Excellent material compatibility, also suitable for flexible endoscopes
- Active against polyomavirus SV 40



- VAH<sup>1</sup> listed; included in the IHO<sup>4</sup> list of virucidal products
- Compatible with all common disinfecting active substances
- Especially user-friendly, pleasant odour
- Does not fixate proteins
- Suitable for stainless steel, anodized aluminium, non-ferrous metals (copper and brass), silicone, synthetic materials, hard plastic, hard rubber, ceramics
- Especially suitable for use in gynaecology

<sup>1</sup> Verbund für Angewandte Hygiene [Association of Applied Hygiene] the Control of Virus Diseases] <sup>3</sup> Carcinogenic, Mutagenic, toxic to Reproduction

<sup>2</sup> Testing guidelines of the Robert Koch Institute and the Deutsche Vereinigung zur Bekämpfung von Viruskrankheiten [German Association for

<sup>4</sup> Industrieverband Hygiene und Oberflächenschutz [German Association for Hygiene and Surface Protection]