



DR. WEIGERT
Systematic Hygiene

Medical | Home Care | Laboratory



neoform[®]
Classic

Amine-based Surface Disinfectant for Surfaces in the Medical, Home Care and Laboratory Sectors.

neoform® Classic –

Amine-based Surface Disinfectant for Surfaces in the Medical, Home Care and Laboratory Sectors.

neoform® Classic¹

For daily combined surface cleaning and disinfection via wiping of the following surfaces:

- Washable surfaces of medical devices and medical equipment as well as other surfaces in hospitals and surgeries
- Washable surfaces in geriatric and care institutions, laboratory areas and public areas
- Bactericidal, yeasticidal and tuberculocidal action, active against enveloped viruses, limited spectrum of virucidal activity, confirmed by experts
- Suitable for materials such as stainless steel, anodised aluminium, plastics (including silicone) and PVC floors

Application:

- Depending on the desired action, dilute a solution with water as per the concentration stated below
- Wipe the surface to be treated with a sufficient quantity of application solution and ensure thorough wetting
- Active substances in 100 g: 9 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine, 0,9 g N-(3-aminopropyl)-N-alkylglycine

Compatibility:

- Application in the **neoform® wipes RTF** wipes dispenser meets the requirements of the current VAH² recommendation
- The instructions for use in the **neoform® wipes RTF** wipes dispenser can be found in the data sheet **neoform® wipes RTF**
- If the dispenser system is closed properly, the soaked neoform wipes can be used for 28 days

Application recommendation (20 °C):

Disinfecting activity
Dirty conditions (non pre-cleaned surfaces)

Bactericidal (VAH ² , EN 13727, EN 16615)	5 ml/l (0.5%), 5 min
Yeasticidal (VAH ² , EN 13624, EN 16615)	5 ml/l (0.5%), 5 min
Tuberculocidal (EN 14348)	15 ml/l (1.5%), 60 min
Virucidal against enveloped viruses (including HIV, HBV, HCV) (EN 14476)	10 ml/l (1.0%), 15 min 5 ml/l (0.5%), 30 min
Limited spectrum virucidal activity (norovirus, adenovirus, including all enveloped viruses) (EN 14476)	20 ml/l (2.0%), 120 min
Rotaviruses (EN 14476)	20 ml/l (2.0%), 90 min
VAH ² list	5 ml/l (0.5%), 5 min



5 l canisters
4470 35



6 x 2 l bottles
4470 42

Our Competence – Your Benefits:

- Short exposure times, broad activity spectrum
- Free from aldehydes and quaternary ammonium compounds
- Good cleaning performance, non protein-fixing, outstanding material compatibility
- User-friendly



- Fast and safe disinfection, thus saving time and costs
- Cleaning results with no stickiness
- Effective, gentle and streak-free surface disinfection
- Pleasant odour, precise and simple dosing

¹ Use disinfectant safely. Always read the label and product information before use.

² Verbund für Angewandte Hygiene [Association of Applied Hygiene] "Recommendation for checking critical points when using wipes dispenser systems in a pre-soak system for surface disinfection" Hyg Med 2012; 37 – 11468