

neodisher® Entschäumer S







Defoamer

Liquid concentrate

Fields of application:

- Preventing foam in special washers for laboratory glassware in medical laboratories (for blood tests and in serology), and laboratories in industry, the food industry and the cosmetics industry, as well as for automated reprocessing of surgical tables
- Also suitable for defoaming in professional warewashers

Performance spectrum:

- Prevents the development of foam during the removal of foam-generating residues such as surfactants, proteins or soaps
- Suitable for use as an additional component with all neodisher detergents for the preliminary or main cleaning
- Based on emulsifying silicone oils, excellent defoaming action even at low temperatures

Special properties:

- Universally applicable
- Highly concentrated and economical
- Odour-neutral

Application and dosage:

neodisher Entschäumer S is used in washer disinfectors for laboratory glassware and decontamination systems for surgical tables as well as in professional dishwashers in addition to neodisher detergents for preliminary cleaning and main cleaning. Depending on the quantity and the kind of foaming soiling neodisher Entschäumer S is dosed at 0.05 – 0.1 ml/l via suitable dosing devices or manually.

Avoid overdosing, otherwise precipitation and deposition of silicone residues on the surfaces of the items being cleaned may result.

neodisher Entschäumer S is dosed either during or immediately after the water intake in washer disinfectors with water exchange. In the case of tank machines a suitable metering pump ensures continuous dosing.

Notes on application:

- For professional use only.
- Do not mix with other products.
- Rinse out dosing system including suction hose with water before changing product
- The neodisher Entschäumer S working solution has to be rinsed off completely (preferably with deionised water).
- Reprocessing of surgical tables:
 - Please observe the reprocessing recommendations of the medical device manufacturers according to the requirements of the DIN EN ISO 17664
 - Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes
- The instructions given by the manufacturer of the washer disinfector and the decontamination system are to be observed
- Storage and transport is only permitted in original packaging.

Technical data:

pH-range	7.0 (concentrate, 20 °C)
Density	approx. 1.0 g/cm³ (20 °C)
Viscosity	approx. 1,000 mPa s (concentrate, 25 °C)





neodisher® Entschäumer S

Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:

< 5 % non-ionic surfactants also preservation agents (benzyloxy)methanol

CE-marks:



neodisher Entschäumer S complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

Storage information:

Always store at a temperature between 5 °C and 30 °C. Usable for 1 year when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol \supseteq .

Hazard and precautionary statements:

neodisher Entschäumer S is not a hazardous product according to the CLP-Regulation (EC) No 1272/2008.

If applied according to the instructions for use the product is safe according to the appropriate guidelines for food processing.


Dispose only when container is empty and closed. For disposal of product residues, refer to Safety Data Sheet.

For further safety information see Safety Data Sheets. These are available at www.drweigert.com under the headline "Service/Downloads".

MB 4301 GB/2-2 Date of issue 03/2021



Tel.: +49-40-7 89 60-0



With the above information, which is appropriate to our current knowledge we describe our product regarding possible safety necessities, but we do not involve any quality description or