




neodisher[®] MA Dental



Alkaline detergent for the automated reprocessing of dental instruments



Powder

Fields of application:

- Automated cleaning of dental instruments and dental utensils such as trays, cups for mouth washing and kidney dishes

Performance spectrum:

- Good cleaning and emulsifying capacity for blood and protein residues
- Protects dental instruments and dental utensils made of stainless steel against corrosion
- Not suitable for the automated cleaning of dental handpieces. For this purpose we recommend the detergent neodisher MediClean Dental

Special properties:

- High material protection
- Particularly suitable for programs with thermal disinfection

Application and dosage:

neodisher MA Dental is used in washer disinfectors for dental instruments. neodisher MA Dental is dosed via the door dosing device for powder and is 3 g/l.

Residues of composite, dental cement and compounds for dental imprints must be removed from the instrument before they set.

For avoiding water stains the use of deionised water in the final rinse is recommended.

The neodisher MA Dental-solution has to be rinsed off completely (preferably with deionised water).

Notes on application:

- For professional use only
- Do not mix with other products

- Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes
- Please observe the reprocessing recommendations of the medical device manufacturers according to the requirements of the DIN EN ISO 17664.
- The instructions given by the manufacturer of the washer disinfectant are to be observed

Technical data:

pH value	11.5 (3 g/l, determined in deionised water, 20 °C)
Bulk density	approx. 1.015 - 1.065 g/l
Titration factor	0.54 (in accordance with the neodisher titration method)

Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:

- < 5 % non-ionic surfactants
- > 30 % phosphates

CE-mark:

neodisher MA Dental complies with European guidelines for medical devices.


If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

Storage information:

Keep the container tightly closed. The product tends to get lumpy when exposed to dampness which can cause a loss of effectiveness. Always store at a temperature between 0 °C and 30 °C. Usable for 2 years when stored as recommended.



neodisher[®] MA Dental

For expiry date refer to the stamp mark on the label behind the hourglass symbol .

Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available e.g. at www.drweigert.com under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

MB 4113/3-2
Date of issue: 02/2021

With the above information, which is appropriate to our current knowledge we describe our product regarding possible safety necessities, but we do not involve any quality description or promise certain properties.
