




# neodisher® MediKlar special



Rinse aid for the automated reprocessing of thermostable and thermolabile instruments



## Liquid concentrate

### Fields of application:

- Rinsing of medical devices such as surgical instruments incl. ophthalmological instruments, anaesthesia equipment, sterilisation containers, and implants as well as baby bottles in washer disinfectors
- Rinsing of bedframes in decontamination systems whenever the bedframes are suitable for automated reprocessing according to the manufacturer's instructions for use

### Performance spectrum:

- Excellent wetting properties, even on hard-to-wet plastic surfaces
- For quick and streak-free drying of the items to be washed
- When neodisher MediKlar special is used as intended, the biocompatibility of the treated instruments is not affected. (Biological evaluation of medical devices according to EN ISO 10993-1)
- Accessories of anaesthesia utensils made of polysulfone (PSU) and polyphenylsulfone (PPSU) e.g. connectors of larynx masks, valves, adapters and parts of instrument container lids can form tension cracks due to the manufacturing process when treated with rinse aids. These must be tested for suitability before use

### Special properties:

- Reduces the drying time considerably
- Neutralizes carried over alkali residues
- Reduces the formation of stains when using softened water in the final rinse

### Application and dosage:

neodisher MediKlar special is used in washer disinfectors and decontamination systems. The

dosage depends among other things on the respective field of application, the items to be washed and the water quality. The following parameters must be observed when using neodisher MediKlar special:

- For the rinsing of surgical instruments, anaesthetic equipment, sterilisation containers, implants and baby bottles: 0.2 - 0.4 ml/l
- For the rinsing of bedframes: 0.2 - 0.8 ml/l

Washer disinfectors, decontamination systems and dosing systems must be suitable for the use of neodisher MediKlar special.

### Notes on application:

- For professional use only.
- Do not mix with other products.
- Rinse out dosing system including suction hose with water before changing product.
- Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes.
- Please observe the reprocessing recommendations of the medical device manufacturers according to the requirements of the DIN EN ISO 17664, the recommendations for reprocessing given by the bedframe producers as well as the recommendations of the German working group "Bedframe and Cart Decontamination Systems" (in the current issue of the AK-BWA brochure "Automated Decontamination").
- The instructions given by the manufacturer of the washer disinfectant are to be observed.



# neodisher<sup>®</sup> MediKlar special

## Technical data:

pH-value	4.0 - 3.5 (0.2 - 0.8 ml/, determined in deionised water, 20 °C)
Viscosity	< 50 mPa s (concentrate, 20 °C)
Density	approx. 1.05 g/cm <sup>3</sup> (20 °C)

## Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:


< 5 % phosphonates, 15 - 30 % non-ionic surfactants, also preservatives (methylchloroisoithiazolinone/ methylisothiazolinone)

## CE-mark:

neodisher MediKlar special complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

## Storage information:

Always store at a temperature between 0 °C and 30 °C. Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

## Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available e.g. at [www.drweigert.com](http://www.drweigert.com) under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

MB 4073/3-2  
Date of issue: 07/2020

With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.