




neodisher®



Active oxygen-based cleaning enhancer for the automated cleaning of instruments using the Miele OXIVARIO process

Liquid concentrate

Fields of application:

Use as a component in combination with alkaline detergents in the "Oxivario" and "Orthovario" processes developed and recommended by Miele

Application and dosage:

Please observe the instructions of the Miele company concerning use and dosage as well as material tolerance.¹⁾

After the whole Miele OXIVARIO-process it is recommended that hinge instruments be treated with neodisher IP Spray.

Notes on application:

- The neodisher OXIVARIO solution has to be rinsed off completely (preferably with deionised water).
- Do not mix with other product concentrates.
- For avoiding water stains the use of deionised water in the final rinse is recommended.
- Use suitable dosing devices.
- For professional use only.
- Rinse out dosing system including hose with water before changing product.
- Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes.
- Please observe the reprocessing recommendations of the medical device manufacturers according to the requirements of the DIN EN ISO 17664.

Technical data:

pH range	5.7 - 5.5 (0.3 - 0.7 %, determined in deionised water, 20 °C)
Density	approx. 1.1 g/cm ³ (20 °C)
Viscosity	< 50 mPas

Ingredients:


Ingredients according to Regulation (EC) No 648/2004 on detergents:
> 30 % oxygen-based bleaching agents

CE-mark:

neodisher OXIVARIO complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

Storage information:

Always store at a temperature between - 20 °C and 25 °C. Keep away from sunlight. Usable for 1 year when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available e.g. at www.drweigert.com under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to Safety Data Sheet.

MB 4004/3-2
Date of issue: 09/2020

¹⁾ The Chemische Fabrik Dr. Weigert GmbH & Co. KG is not responsible for any legal obligations (guarantees) concerning possible changes of medical devices reprocessed with these Miele processes.

⁹⁾ With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.