

# Rinse aid for the automated reprocessing of animal cages





## Fields of application:

Rinsing of animal cages and accessories as well as cage, storage and transport racks during automated reprocessing

#### Performance Spectrum:

- Improves and accelerates drying considerably, particularly of plastic materials
- Suitable for materials usually used for animal cages, e.g. polysulfones, polycarbonates, polyphenylsulfones and polyetherimides

#### **Special Properties:**

- Free of surfactants, acidic formulation neutralises alkaline residues
- Excellent material compatibility, prevents the formation of stress cracks during autoclaving
- Toxicologically safe

## Application and dosage:

neodisher PolyKlar is used in specially designed washers for animal cages and accessories as well as cage, storage and transport racks. The dosing amount depends among other things on the material of the items to be washed and is 0.5 - 2.0 ml/l.

Use suitable dosing devices.

For the automated cleaning of animal cages and accessories we recommend using neodisher LaboClean FLA in the alkaline process and neodisher N in the acidic process.

#### General instructions on use:

- Do not mix with other products.
- Rinse out dosing system including hose with water before changing product.
- For professional use only.
- The instructions given by the manufacturer of the washer are to be observed.
- Please observe the recommendations on reprocessing given by the manufacturers of the animal cages and the items to be washed as well as the recommendations of the working group for cage processing (AK KAB)

in the current issue of the AK KAB brochure "Cage Processing in Animal Facilities properly done".

## Technical data:

pH-range	3.6 - 3.2 (0.5 - 2 ml/l, determined
	in deionised water, 20 °C)
Density	approx.1.0 g/cm³ (20 °C)
Viscosity	< 10 mPa s (concentrate, 20 °C)

### Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents: < 5 % polycarboxylates, also: preservatives (chloromethylisothiazolinones/methylisothiazolinones)

#### Storage information:

Always store at a temperature between 0 °C and 30 °C. Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol  $\underline{\square}$ .

#### Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available e.g. at www.drweigert.com under the headline "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to Safety Data Sheet.

MB 4019/3-2 Revision Date 09/2020

With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.

