

# neodisher<sup>®</sup> SBN

Liquid acidic rinsing agent for use in bedpans washers



<b>Main fields of application:</b>	Used as a rinsing agent and complexing agent for the water hardness in bedpan washers with steam disinfection for the preparation of bedpans, urine bottles and stainless steel bowls and pans in hospitals and care homes.
<b>Characteristics:</b>	<p>neodisher SBN is a slightly acidic rinse aid with special softening components in a low-foam formulation. Because of its good lime bonding capacity and the wetting properties the bedpans are rinsed clean without spotting and without lime deposition. When used in accordance with the instructions there is no lime deposition in the steam generator or in the pipes.</p> <p>neodisher SBN can be used in water of any hardness. Rinse out dosing system including hose with water before changing product.</p>
<b>Application and dosage:</b>	<p>0.5 - 1.5 ml/l depending on water hardness through built-in metering device. Note acid resistance of metering systems and steam producers.</p> <p>With firmly adherent contamination the additional use of neodisher SBR, a special cleaning agent for automatic washing machines for bedpans is recommended.</p> <p>Do not mix with other products. Rinse out dosing system including suction hose with water before changing product.</p>
<b>Technical data:</b>	<p>Density (20 °C): 1.1 g/cm<sup>3</sup> pH-range (determined in deionised water, 20 °C) 0.5 - 1.5 ml/l: 3.7 - 3.0 Viscosity (concentrate, 20 °C): &lt; 50 mPas</p>
<b>Ingredients:</b>	<p>Ingredients according to Regulation (EC) No 648/2004 on detergents: 15 - 30 % non-ionic surfactants also perfumes (linalool), preservation agents (methylchloroisothiazolinone/methylisothiazolinone)</p>
<b>CE-marks:</b>	neodisher SBN conforms to the requirements of the legislation for medical devices according to the European directive 93/42/EC.
<b>Storage information:</b>	The product is sensitive to frost. Usable for 3 years when stored as recommended.
<b>Hazard and precautionary statements:</b>	<p>For safety information see EC safety data sheets. These are available at <a href="http://www.drweigert.com">www.drweigert.com</a> under the category "Download".</p> <p>Dispose only when container is empty and closed. For disposal of product residues, refer to the Material Safety Data Sheet.</p>

MB 4042 GB/2-4  
05/15



Chemische Fabrik Dr. Weigert GmbH & Co. KG  
Mühlenhagen 85, D - 20539 Hamburg

Telefon: (040) 789 60 - 0  
Telefax: (040) 789 60 - 120

E-Mail: [info@drweigert.de](mailto:info@drweigert.de)  
Internet: [www.drweigert.de](http://www.drweigert.de)

Die Angaben dieses Merkblattes beruhen auf unseren derzeitigen Kenntnissen und Erfahrungen. Sie befreien den Verwender nicht von eigenen Prüfungen und Versuchen. Eine rechtlich verbindliche Zusicherung bestimmter Eigenschaften kann hieraus nicht abgeleitet werden.

With the above information, which is appropriate to our current knowledge we describe our product regarding possible safety necessities, but we do not involve any quality description or promise certain properties.