




# neodisher<sup>®</sup> Septo PreDis ZP



## Detergent and pre-disinfectant for thermostable and thermolabile instruments



### Liquid concentrate

#### Fields of application:

- Disinfecting pre-cleaning of thermostable and thermolabile instruments, dental instruments and flexible endoscopes e.g. prior to automated reprocessing for optimum personnel safety
- Wet disposal of surgical instruments and dental instruments
- Also suitable for the manual cleaning and disinfection of thermostable and thermolabile instruments including flexible endoscopes

#### Performance spectrum:

- Bactericidal, yeasticidal and active against enveloped viruses (incl. HBV, HCV and HIV)
- The disinfecting activity has been proved in accordance with European methods and VAH<sup>1</sup> methods
- The activity against enveloped viruses has been tested and confirmed in accordance with European methods
- Suitable for stainless steel, anodised aluminium, silicone, plastics
- VAH<sup>1</sup>-listed, included in the IHO<sup>2</sup> list of virucidal products

#### Special properties:

- Excellent cleaning performance, does not fixate proteins
- Excellent material compatibility
- Free of aldehydes

#### Application and dosage:

neodisher Septo PreDis ZP can be used in immersion baths or ultrasonic baths. Depending on the desired activity a solution is prepared according to the dosing recommendation given below. Place the instruments or endoscopes in this solution and observe the below-mentioned contact times.

To prepare neodisher Septo PreDis ZP solutions it is recommended to use softened water or soft

water with less than 3 °d (German hardness). The water hardness should not exceed 20 °d (German hardness).

Application recommendation (20 °C, dirty conditions)	
bactericidal (EN 13727 and EN 14561) yeasticidal (EN 13624 and EN 14562)	5 ml/l (0.5 %), 15 min
bactericidal/yeasticidal (VAH <sup>1</sup> )	5 ml/l (0.5 %), 15 min
active against enveloped viruses (incl. HIV, HBV, HCV) (EN14476, EN 17111)	5 ml/l (0.5 %), 15 min

The neodisher Septo PreDis ZP solution is to be rinsed off with water (preferably deionised).

The shelf life of the unused application solution is 7 days.

We recommend renewing the used solution daily; in the event of a higher dirt load more frequently.

#### General notes on application:

- For professional use only.
- For controlled dosage the use of manual dosing aids, e.g. dosing devices for disinfectants is recommended. Please contact us.
- It is generally recommended to wear gloves when working with disinfectants.
- Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes.
- Please observe the reprocessing recommendations of the medical device manufacturers according to the requirements of the DIN EN ISO 17664
- Do not mix with other products



# neodisher<sup>®</sup> Septo PreDis ZP

## Expert reports:

The disinfecting activity has been tested and confirmed by certification. We will be pleased to provide certificates on request.

## Technical data:

pH-value	9.1 (5 ml/l, determined in deionised water, 20 °C)
Viscosity	< 50 mPa s (concentrate, 20 °C)
Density	approx. 1.0 g/cm <sup>3</sup> (20 °C)

Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

MB 4208/3-5  
04/2020

1 Verbund für Angewandte Hygiene [Association of Applied Hygiene]  
2 Industrieverband Hygiene und Oberflächenschutz [Industry Association Surface Protection and Hygiene, Germany]

## Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:

5 - 15 % non-ionic surfactants

Active substances in 100 g:

10.0 g didecyldimethylammoniumchloride,  
2.0 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine

CE-mark: 

neodisher Septo PreDis ZP complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

## Storage information:

Always store at a temperature between 0 °C and 30 °C.

Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

## Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available at [www.drweigert.com](http://www.drweigert.com) under the category "Service/Downloads".



With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties..