




neodisher® SystemAct



Enzymatic cleaning component
for processing thermostable
and thermolabile instruments

Highly concentrated liquid

Fields of application:

Automated cleaning of surgical instruments, incl. MIS instruments and micro instruments, anaesthesia utensils, containers and other medical utensils

Performance spectrum:

Enzymatic cleaning component which is dosed together with the mildly alkaline cleaning component neodisher SystemClean into the cleaning step of the washer disinfectant. The resulting enzymatic alkaline solution has the following properties:

- Reliably removes residues of dried and denatured blood; at the same time exhibiting excellent material compatibility
- Fulfills the current recommendations of the Robert Koch Institute (RKI) for processing medical devices for minimizing the risk of transmitting the new variant of the Creutzfeldt Jakob Disease (vCJD)
- Supports the removal of biofilms
- Suitable for instruments and utensils made of stainless steel, instrument steel, anodised aluminium, optics, commonly used synthetic materials and materials of anaesthesia utensils
- Excellent compatibility with anodised aluminium when dosed according to individual requirements

Special properties:

- Component concentrations can be individually selected to meet the requirements, depending on soiling
- Highly concentrated liquid: maximum efficiency
- Excellent material compatibility
- Low-foam formula
- First-class cleaning performance due to the individually adaptable combination of alkali and enzyme components

- No neutralisation step necessary with automated processing, therefore short program cycles
- Fitted with RFID transponder for automatic product identification with the weigomatic system ALPHA component dosing system to safely exclude product mix-ups

Application and dosage:

The cleaning components neodisher SystemClean and neodisher SystemAct are used in washer disinfectants. The dosing amount depends among other things on the area of application, the items to be cleaned, the wash mechanics and the degree of soiling of the instruments. The following parameters are recommended when using neodisher SystemClean and neodisher SystemAct:

| | |
|---|---|
| Automated cleaning of thermostable and thermolabile instruments | 0.5 - 1 ml/l (0.05 - 0.1 %) neodisher SystemAct |
| | 0.5 - 1 ml/l (0.05 - 0.1 %) neodisher SystemClean |
| | e.g. 30 - 60 °C, 10 min |

The RFID-secured, mix-up-free provision of the components through storage containers is carried out via the weigomatic system ALPHA which has been especially developed for dosing highly concentrated liquids.

Suitable dosing devices are to be used in the washer disinfectants.

The use of deionised water is recommended in the cleaning step and the final rinse. When using deionised water in the final rinse water stains are avoided and anodised aluminium is protected at the same time.



DR. WEIGERT

Chemische Fabrik Dr. Weigert GmbH & Co. KG
Mühlenhagen 85 • D-20539 Hamburg

Tel.: +49-40-7 89 60-0
Fax: +49-40-7 89 60-120

info@drweigert.de
www.drweigert.com



neodisher[®] SystemAct

The neutralisation step required when using classic alkaline detergents can be omitted. When processing ophthalmological instruments an additional intermediate rinse is recommended prior to the final rinse.

Notes on application:

- For professional use only.
- Do not mix components as concentrates.
- Do not mix with other products.
- Rinse out dosing system including suction hose with water before changing product.
- Processing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes.
- The solution has to be rinsed off completely (preferably with deionised water).
- Please observe the medical device manufacturers recommendations for use according to the requirements of DIN EN ISO 17664.
- The instructions given by the manufacturer of the washer disinfectant are to be observed.

Technical data:

| | |
|-----------|--|
| pH range | 6.5 – 6.7 (0.5 – 1.0 ml/l, determined in deionised water, 20 °C) |
| Density | approx. 1.0 g/cm ³ (20 °C) |
| Viscosity | < 10 mPa s (concentrate, 20 °C) |

Ingredients:


Ingredients according to Regulation (EC) No 648/2004 on detergents:
< 5 % non-ionic surfactants, polycarboxylates
also: enzymes, preservatives (2-phenoxyethanol)

CE-mark:

neodisher SystemAct complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

Storage information:

Store in a cool but frost-free place. Always store at a temperature between 0 °C and 25 °C. Keep away from sunlight. Usable for 2 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available at www.drweigert.com under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

MB 4023/2-8

Date of Issue 09/2025

With the above information, which is appropriate to our current knowledge we describe our product regarding possible safety necessities, but we do not involve any quality description or promise certain properties.