

## neodisher<sup>®</sup> TN

# Neutral special rinse aid for difficult-to-wet materials

#### Liquid concentrate

#### Fields of application:

- Rinsing dishes, glasses, cutlery, trays and utensils in professional dishwashers as well as containers and trolleys in decontamination units in all food processing areas, such as professional kitchens, butcher's shops and bakeries
- Rinsing difficult-to-wet materials such as bed frames, bedside cupboards and trolleys in decontamination units in clinics and in the care sector

#### Performance spectrum:

- Achieves a brilliant rinsing result due to excellent wetting of all materials
- Good material compatibility, specially matched to the various materials of bedframes and bedside cupboards
- Especially suitable for use with soft water
- When neodisher TN is used as intended, the biocompatibility of the treated medical devices is not affected. (Biological evaluation of medical devices according to EN ISO 10993-1)

#### Special properties:

- Achieves excellent wetting and stain and streak-free drying; also of difficult-to-wet parts made of synthetic materials
- Reduced consumption
- Even and residue-free drying

#### Application and dosage:

neodisher TN can be used in professional dishwashers and decontamination units. neodisher TN is dosed into the final rinse water via an automatic dosing device. Dosage depends on water quality, shape and material of the items to be washed and the temperature of the rinse water. The optimum dosage of neodisher TN should be adjusted when commissioning the washer.





Depending on water quality,<br/>items to be washed and<br/>temperature of the rinse water0.1 - 0.8 ml/l

When working with hard or salty water a special water treatment is recommended to achieve an optimal rinsing result. In order to avoid water stains we recommend using deionised water in the final rinse step.

#### Notes on application:

- For professional use only.
- Do not mix with other products.
- Rinse out dosing system including suction hose with water before changing product.
- The instructions given by the manufacturer of the dishwasher are to be observed.
- Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes.
- Please observe the reprocessing recommendations of the medical device manufacturers according to the requirements of the DIN EN ISO 17664, the reprocessing recommendations of the manufacturer of the bedframes and the recommendations of the Working Group Bedframe and Trolley Decontamination Units [Arbeitskreis Bettgestell- und Wagendekontaminationsanlagen (AK-BWA)] in the current issue of the AK-BWA brochure "Automated Decontamination" ["Maschinelle Dekontamination"].
- Storage and transport is only permitted in original packaging.





## neodisher<sup>®</sup> TN

#### Technical data:

pH range	6.7 – 7.0 (0.1 – 0.8 ml/l, determined
	in deionised water, 20 °C)
Viscosity	< 50 mPa s
	(concentrate, 20 °C)
Density	approx. 1.0 g/cm <sup>3</sup> (20 °C)

#### Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:

< 5 % phosphonates, polycarboxylates 15 - 30 % non-ionic surfactants also preservatives: Reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one [EC No. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC No. 220-239-6] (3:1)

### CE mark: CE MD

neodisher TN complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

#### Storage information:

Always store at a temperature between 0 °C and 30 °C. Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol  $\blacksquare$ .

#### Hazard and precautionary statements:

neodisher TN is not a hazardous product according to the CLP Regulation (EC) No 1272/2008.

If applied according to the instructions for use the product is safe according to the appropriate guidelines for food processing.

Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

For further safety information see Safety Data Sheets. These are available at www.drweigert.com under the category "Service/Downloads".

> MB 3115/3-4 Date of issue: 06/2021

With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.

