



neiform MED AF



Disinfecting cleaner for surfaces

Liquid concentrate

Fields of application:

- Cleaning and disinfecting surfaces of medical devices, medical equipment and other surfaces in hospitals, doctors' practices, old people's homes and laboratories
- Also suitable for use with the refillable wiping system Kimtech* Wettask* DS of Kimberly-Clark

Performance spectrum:

- Bactericidal and yeasticidal activity according to DGHM¹-methods. Also active against enveloped viruses according to RKI/DVV².
- The application in the bucket wipes system Kimtech* Wettask* DS of Kimberly-Clark fulfils the requirements of the current VAH recommendation³
- VAH⁴-listed
- included in the IHO⁵ list of virucidal products
- Free of aldehydes

Special Properties:

- Wide activity spectrum, even at low concentrations resp. short action times
- Excellent cleaning performance
- Excellent material compatibility.

Application and dosage:

neiform MED AF is used for the combined cleaning and disinfection of surfaces at room temperature via wiping. For achieving the desired activity dilute neiform MED AF according to the below-mentioned dosing recommendations. Wipe the surfaces with the working solution. Adhere to the below-mentioned contact times.

neiform MED AF can also be used with the refillable wiping system Kimtech* Wettask* DS of Kimberly-Clark: Place the roll of wipes into the wipes dispenser system. Fill from the prepared working solution (with a concentration of at least 0.5%) 2.0 l into the bucket or 0.75 l into the box.

Pull the first wipe (in the middle of the wipes roll) through the opening of the dispenser and tightly

close the wipes dispenser system. After 60 min the wipes are completely soaked and ready for use. Simply pull out the soaked wipes through the opening of the dispenser in the required amount. After removal tightly close the dispenser system. The soaked wipes can be used for up to 28 days when the wiping system is kept tightly closed.

The VAH recommends³ thorough disinfecting cleaning and drying after depletion and before replenishing the dispenser system. The wipes system Kimtech* Wettask* DS of Kimberly-Clark can be reprocessed manually as well as thermally (thermostable till 100 °C) or chemothermally in a washer-disinfector.

Disinfecting performance (20°C) clean and dirty conditions (pre-cleaned and non-pre-cleaned surfaces)	
bactericidal / yeasticidal (DGMH)	20 ml/l (2.0 %), 5 min 10 ml/l (1.0 %), 15 min 5 ml/l (0.5 %), 30 min 2.5 ml/l (0.25 %), 60 min
active against enveloped viruses(RKI/DVV) (incl. HIV, HBV, HCV)	20 ml/l (2.0 %), 5 min 10 ml/l (1.0 %), 15 min 5 ml/l (0.5 %), 30 min

After cleaning and disinfecting surfaces must be rinsed with drinking water in order to avoid residues. This relates in particular to powder-coated and stainless steel surfaces as these may cause changes in the material.

General instructions on use:

- For professional use only!
- For the economical and controlled dosage, the use of manual dosing aids such as mixing equipment or dosing caps is recommended. Please contact us.
- It is generally recommended to wear gloves when working with disinfectants.
- After cleaning and disinfecting measures surfaces must be rinsed with drinking water in order to avoid residues.
- Reprocessing should comply with all ordinances pursuant to the Medical Device Directive and



neiform MED AF

should be performed with appropriate validated processes.

- Do not mix with other products.
- *Use disinfectants safely. Always read the label and product information before use.*

Expert Reports / listings:

The disinfecting activity as well as the service life and the compatibility in the wipes dispenser system Kimtech* Wettask* DS of Kimberly-Clark have been tested and confirmed by certification in accordance with the VAH recommendation³. We will be pleased to provide certificates on request.

Technical data:

pH-value	approx. 6 (concentrate, 20 °C)
Density	1.0 g/cm ³ (20 °C)
Viscosity	< 10 mPa s

Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:

5 - 15 % non-ionic surfactants
also disinfectants, perfumes


active substances in 100 g:
17.5 g didecyldimethylammonium chloride
4.0 g poly(hexamethylene biguanide)

CE-marks:  0297

neiform MED AF conforms to the European directive 93/42/EC, Annex I, concerning medical devices.

Storage information:

Always store at a temperature between 0 °C and 25 °C. Keep container tightly closed. Usable for 4 years when stored as recommended.

For expiry date refer to the stamp mark on the label behind the hourglass symbol .

Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available at www.drweigert.com under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

MB 4056/3-6
Revision date 06/2016

- 1 according to the guidelines of the Deutschen Gesellschaft für Hygiene und Mikrobiologie [DGHM, German Society of Hygiene and Microbiology]
- 2 according to the guidelines of the Robert Koch-Institute (RKI) and the Deutschen Vereinigung zur Bekämpfung von Viruskrankheiten [DVV, German Association for the control of Virus Diseases]
- 3 recommendation for checking critical points when using wipes dispenser systems in a pre-soak system for surface disinfection (Hyg Med 2012; 37-11)
- 4 Verbund für Angewandte Hygiene [Association of Applied Hygiene]
- 5 Industrieverband Hygiene und Oberflächenschutz [Industry Association Surface Protection and Hygiene, Germany]



With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.