




neoform MED FF



Disinfectant cleaner for surfaces

Liquid concentrate

Fields of application:

- Cleaning and disinfecting surfaces of medical devices, medical equipment and other surfaces in hospitals, in doctors' practices, old people's homes as well as laboratories.
- Also suitable for use with the Dr. Weigert fleece wipe dispensing system neoform wipes RTF

Performance spectrum:

- Certified bactericidal, yeasticidal and tuberculocidal activity according to DGHM¹-methods, virucidal activity according to RKI/DVV²-methods as well as active against Clostridium difficile according to EN method
- VAH³-listed
- Included in the IHO⁴ disinfectant list
- The application in the Dr. Weigert fleece wipe dispensing system neoform wipes RTF fulfils the requirements of the current VAH recommendation⁵

Special properties:

- Wide activity spectrum and good cleaning performance
- Also suitable for disinfection in the event of an outbreak of C. difficile
- Free of formaldehyde

Application and dosage:

neoform MED FF is used for combined cleaning and disinfection of surfaces at room temperature via wiping. For achieving the desired activity, dilute neoform MED FF according to the below-mentioned dosing recommendations.

Wipe the surfaces to be cleaned and disinfected with the solution and adhere to the below-mentioned contact times.

For instructions on how to use the product with the neoform wipes RTF fleece wipe dispensing system please refer to the neoform wipes RTF data sheet. The soaked neoform wipes can be used for up to 28 days when the wipes system is kept tightly closed.

Disinfecting performance (20 °C)		
	clean conditions (pre-cleaned surfaces)	dirty conditions (non-pre-cleaned surfaces)
bactericidal / yeasticidal (DGHM)	10 ml/l (1.0%), 5 min	15 ml/l (1.5%), 5 min
	5 ml/l (0.5%), 15 min	10 ml/l (1.0%), 15 min
	2.5 ml/l (0.25%), 60 min	5 ml/l (0.5%), 30 min 2.5 ml/l (0.25%), 60 min
tuberculocidal (DGHM)	40 ml/l (4.0 %), 30 min	40 ml/l (4.0 %), 30 min
	10 ml/l (1.0 %), 60 min	10 ml/l (1.0 %), 60 min
active against enveloped viruses (RKI/DVV) (incl. HIV, HBV, HCV)	5 ml/l (0.5 %), 5 min	-
virucidal (RKI/DVV)	20 ml/l (2.0 %), 60 min	-
active against C. difficile (EN 13704)	30 ml/l (3.0 %), 30 min	-

When used on surfaces of medical devices a rinsing step or a wiping step with drinking water is necessary after the contact time has expired. This relates in particular to powder-coated and stainless steel surfaces as these may cause surface changes in the material.

General instructions on use:

- For professional use only!
- For the economical and controlled dosage the use of manual dosing aids such as mixing



neoform MED FF

equipment or dosing cups is recommended.
Please contact us.

- It is generally recommended to wear gloves when working with disinfectants.
- Reprocessing should comply with all ordinances pursuant to the Medical Device Directive and should be performed with appropriate validated processes.
- Do not mix with other products.
- *Use disinfectants safely. Always read the label and product information before use.*

Expert reports:

The disinfecting activity has been confirmed by certification. We will be pleased to provide certificates on request.

Technical data:

pH-value	4.5 (concentrate, 20 °C)
Viscosity	< 10 mPa s
Density	approx. 1.0 g/cm ³ (20 °C)

Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:
< 5 % non-ionic surfactants
also disinfectants

Active substances in 100 g:


5.0 g benzalkonium chloride
3.0 g didecyldimethylammonium chloride
8.0 g glutaral

CE-mark:  0297

neoform MED FF complies with European guidelines for medical devices.

Storage information:

Always store at a temperature between 0 °C and 25 °C. Keep container tightly closed.
Usable for 4 years when stored as recommended.

For expiry date refer to the stamp mark on the label behind the hourglass symbol .

Hazard and precautionary statements:

For safety information see EC safety data sheets. These are available at www.drweigert.com under the category "Service/Downloads".

Dispose only when container is empty and closed. For the disposal of product residues, refer to Safety Data Sheet.

MB 4057/3-2
Date of Issue: 11/2019

- 1 according to the guidelines of the Deutschen Gesellschaft für Hygiene und Mikrobiologie [DGHM, German Society of Hygiene and Microbiology]
- 2 according to the guidelines of the Robert Koch-Institute (RKI) and the Deutschen Vereinigung zur Bekämpfung von Viruskrankheiten [DVV, German Association for the control of Virus Diseases]
- 3 Verbund für Angewandte Hygiene [Association of Applied Hygiene]
- 4 Industrieverband Hygiene und Oberflächenschutz [Industry Association Surface Protection and Hygiene, Germany]
- 5 recommendation for checking critical points when using wipes dispenser systems in a pre-soak system for surface disinfection (Hyg Med 2012; 37-11)



With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.