



Disinfectant cleaner for surfaces

Liquid concentrate

Fields of application:

- Cleaning and disinfecting surfaces of medical devices, medical equipment and other surfaces in hospitals, in doctors' practices, old people's homes as well as laboratories.
- Also suitable for use with the Dr. Weigert fleece wipe dispensing system neiform wipes RTF

Performance spectrum:

- Certified bactericidal, yeasticidal and tuberculocidal activity according to VAH¹-methods, virucidal activity according to RKI/DVV²-methods as well as active against Clostridium difficile according to EN method
- VAH¹-listed
- The application in the Dr. Weigert neiform wipes RTF fleece wipe dispensing system fulfils the requirements of the current VAH recommendation³

Special properties:

- Broad activity spectrum and good cleaning performance
- Also suitable for disinfection in the event of an outbreak of C. difficile
- Free of formaldehyde

Application and Dosage:

neiform MED FF is used for combined cleaning and disinfection of surfaces at room temperature via wiping. For achieving the desired activity, dilute neiform MED FF according to the below-mentioned dosing recommendations.

Wipe the surfaces to be cleaned and disinfected with the solution and adhere to the below-mentioned contact times.

For instructions on how to use the product with the neiform wipes RTF fleece wipe dispensing system, please refer to the neiform wipes RTF data sheet. The soaked neiform wipes can be used for up to 28 days when the wipes system is kept tightly closed.

Disinfecting performance (20 °C)		
	clean conditions (pre-cleaned surfaces)	dirty conditions (non-pre-cleaned surfaces)
bactericidal / yeasticidal (VAH)	10 ml/l (1.0%), 5 min	15 ml/l (1.5%), 5 min
	5 ml/l (0.5%), 15 min	10 ml/l (1.0%), 15 min
	2.5 ml/l (0.25%), 60 min	5 ml/l (0.5%), 30 min 2.5 ml/l (0.25%), 60 min
tuberculocidal (VAH)	40 ml/l (4.0 %), 30 min	40 ml/l (4.0 %), 30 min
	10 ml/l (1.0 %), 60 min	10 ml/l (1.0 %), 60 min
active against enveloped viruses (RKI/DVV) (incl. HIV, HBV, HCV)	5 ml/l (0.5 %), 5 min	-
virucidal (RKI/DVV)	20 ml/l (2.0 %), 60 min	-
active against C. difficile (EN 13704)	30 ml/l (3.0 %), 30 min	-

When used on surfaces of medical devices a rinsing step or a wiping step with drinking water is necessary after the contact time has expired. This relates in particular to powder-

With the above information, which is appropriate to our current knowledge we describe our product regarding possible safety necessities, but we do not involve any quality description or



coated and stainless steel surfaces as these may cause surface changes in the material

General notes on application:

- For professional use only!
- For the economical and controlled dosage the use of manual dosing aids such as mixing equipment or dosing cups is recommended. Please contact us.
- It is generally recommended to wear gloves when working with disinfectants.
- Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes.
- Do not mix with other products.
- *Use disinfectants safely. Always read the label and product information before use.*

Expert reports:

The disinfecting activity has been confirmed by certification. We will be pleased to provide certificates on request.

Technical data:

pH-value	4.5 (concentrate, 20 °C)
Density	approx. 1.0 g/cm ³ (20 °C)
Viscosity	< 10 mPa s

MB 4057/3-3
Stand: 07/22

Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents: < 5 % non-ionic surfactants also disinfectants

Active substances in 100 g:
5.0 g benzalkonium chloride
3.0 g didecyldimethylammonium chloride
8.0 g glutaral

CE-Kennzeichnung:

neoform MED FF complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

Storage information:

Always store at a temperature between 0 °C and 25 °C. Keep container tightly closed.

Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol

Hazard and precautionary statements:

For safety information see EC safety data sheets. These are available at www.drweigert.com under the category "Service/Downloads".

Dispose only when container is empty and closed. For the disposal of product residues, refer to Safety Data Sheet.

- 1 Verbund für Angewandte Hygiene [Association of Applied Hygiene]
- 2 according to the guidelines of the Robert Koch-Institute (RKI) and the Deutschen Vereinigung zur Bekämpfung von Viruskrankheiten [DVV, German Association for the control of Virus Diseases]
- 3 recommendation for checking critical points when using wipes dispenser systems in a pre-soak system for surface disinfection (Hyg Med 2012; 37-11)