



# neoform MED rapid



## Rapid disinfectant for medical devices

Ready-for-use solution

### Fields of application:

- Specific and rapid disinfection of surfaces of medical devices such as surgical tables in hospitals as well as in doctors' and dentists' practices

### Performance spectrum:

- Bactericidal, yeasticidal, fungicidal and tuberculocidal activity confirmed in accordance with VAH<sup>1</sup> and EN methods<sup>2</sup> as well as activity against enveloped viruses confirmed according to the test guidelines of the RKI/DVV<sup>3</sup> by certification
- Also limited spectrum virucidal activity (active against noro and adenovirus) and active against rotavirus
- The surfaces to be disinfected must be compatible with alcoholic disinfectants
- Not suitable for disinfecting acrylic glass (perspex)
- VAH<sup>1</sup>-listed
- Included in the IHO<sup>4</sup> list of disinfectants

### Special properties:

- Wide activity spectrum, even with short contact times
- Rapid and stain-free drying, no rinsing with fresh water necessary
- Free of aldehydes, colourants and perfumes

### Application and dosage:

Apply neoform MED rapid undiluted on clean and dry surfaces and disinfect by wiping. Surfaces must be completely wetted with neoform MED rapid. To achieve the desired activity, the below-mentioned contact times must be strictly adhered to.

	Clean conditions, 20 °C (cleaned surfaces)
bactericidal (VAH <sup>1</sup> , EN 13727, EN 16615)	undiluted, 1 min
yeasticidal (VAH <sup>1</sup> , EN 13624, EN 16615)	undiluted, 1 min



fungicidal (EN 13624)	undiluted, 20 min
tuberculocidal (EN 14348)	undiluted, 30 sec
Active against enveloped viruses (incl. HIV, HBV, HCV) (EN 14476, RKI/DVV <sup>3</sup> )	undiluted, 15 sec
Limited virucidal activity (adenovirus, norovirus, incl. all enveloped viruses) (EN 14476)	undiluted, 2 min
Active against norovirus (MNV <sup>5</sup> ) (EN 14476)	undiluted, 2 min
Active against adenovirus (EN 14476)	undiluted, 1 min
Active against rotavirus (EN 14476)	undiluted, 1 min
VAH <sup>1</sup> list	undiluted, 5 min

Due to its alcohol content, neoform MED rapid must not be applied in too large amounts (apply max. 50 ml per m<sup>2</sup> on the surfaces to be treated and max. 100 ml per m<sup>2</sup> base area). Please thoroughly read and observe the hazard and precautionary statements.

### General instructions for use:

- For professional use only.
- For economical and controlled dosage the use of manual dosing aids is recommended. Please contact us.
- It is generally recommended to wear gloves when working with disinfectants.
- Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes.
- Please observe the reprocessing recommendations given by the manufacturer of the medical devices and surfaces. Do not mix with other products.

### Expert reports:

The disinfecting activity has been confirmed by certification. We will be pleased to provide certificates on request.



# neoform MED rapid

## Technical data:

Viscosity	< 10 mPa s (undiluted, 20 °C)
Density	approx. 0.9 g/cm <sup>3</sup> (20 °C)
Flash point	27 °C according to DIN 22791

## Ingredients:


Active substances in 100 g:  
35.0 g propan-1-ol,  
25.0 g ethanol



neoform MED rapid complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

## Storage information:

Always store at a temperature between -15 °C and 25 °C. Keep container tightly closed. Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

## Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available at [www.drweigert.com](http://www.drweigert.com) under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to Safety Data Sheet.

MB 4070/3-5 Date of Issue: 11/2020

- 1 Verbund für Angewandte Hygiene [Association of Applied Hygiene, Germany]
- 2 European test methods
- 3 according to the guidelines of the Robert Koch-Institute (RKI) and the Deutschen Vereinigung zur Bekämpfung von Viruskrankheiten (DVV) [German Association for the Control of Virus Diseases]
- 4 Industrieverband Hygiene und Oberflächenschutz [Industry Association Surface Protection and Hygiene, Germany]
- 5 Murine norovirus

With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.