




DR. WEIGERT
Systematic Hygiene

Medical



neodisher®
**Septo
PreClean**

Disinfecting cleaner for the manual reprocessing of thermostable and thermolabile instruments.

neodisher® Septo PreClean – efficient cleaning with increased safety.

neodisher® Septo PreClean

Cleaning and disinfecting agent for optimum personnel safety, with excellent material compatibility

- For disinfecting pre-cleaning of thermostable and thermolabile instruments incl. flexible endoscopes and dental instruments
- For wet disposal of surgical instruments and dental instruments over longer periods of time, e.g. overnight or over the weekend
- Also suitable for the manual cleaning and disinfection of thermostable and thermolabile instruments, incl. flexible endoscopes



5 L canister
4503 35



10x 1 L bottle
4503 45

Application and dosage:

- **neodisher® Septo PreClean** is used in immersion and ultrasonic baths. Depending on the desired activity the working solution must be prepared according to the dosing recommendations.
Place the instruments resp. endoscopes into the working solution and adhere to contact times
- **neodisher® Septo PreClean** is active against bacteria, yeasts and enveloped viruses (incl. HIV, HBV, HCV) when used with the below-mentioned parameters

Application recommendation (20 °C):

Disinfecting pre-cleaning	15 ml/l (1.5 %), 15 min
Wet disposal	10 ml/l (1.0 %), at least 30 min
Disinfection	15 ml/l (1.5 %), 15 min



Our Competence – Your Benefit:

- excellent cleaning performance with disinfecting activity
- certified activity: bactericidal/yeasticidal according to DGHM¹ methods and European standards², also active against enveloped viruses according to RKI/DVV³
- excellent material compatibility; free of aldehydes and quaternary ammonium compounds



- effectively prevents fixing of adhesive protein residues, for an optimized pre-treatment
- improved safety for the personnel handling the instruments, increases hygiene safety associated with the reprocessing, VAH⁴-listed and included in the IHO⁵ list of virucidal products
- universal application: suitable for stainless steel, anodized aluminium, non-ferrous metals (copper and brass), silicone and plastics; very user-friendly

¹ Deutsche Gesellschaft für Mikrobiologie und Hygiene (German Society of Hygiene and Microbiology) ² European Standards EN13727, EN13624, EN14561, EN14562 ³ Robert Koch Institute, Germany/Deutsche Vereinigung zur Bekämpfung der Viruskrankheiten (German Association for the Control of Virus Diseases) ⁴ Verbund für Angewandte Hygiene (Association of Applied Hygiene) ⁵ Industrieverband Hygiene und Oberflächenschutz (Industry Association Surface Protection and Hygiene, Germany)

For further information on products and application please refer to the product data sheet.



Your expert for application queries: Ina Haacke · ina.haacke@drweigert.de