



Always on Hand
to Help with Hygiene –
**in Europe and Around
the World.**



neodisher® –
for Maximum
**Hygiene
Safety.**

Chemische Fabrik
Dr. Weigert GmbH & Co. KG
Mühlenhagen 85
20539 Hamburg
Tel.: +49-40-7 89 60-0
info@drweigert.de

Dr. Weigert
Handelsgesellschaft m.b.H.
Wienerbergstraße 11/12 a
1100 Wien
Tel.: +43-1-99460 6275
info@drweigert.at

Dr. Weigert (Schweiz) AG
Baarerstrasse 14
6300 Zug
Tel.: +41-41 561 32 90
info@drweigert.ch

Dr. Weigert Nederland BV
Balkendwarsweg 3
9405 PT Assen
Tel.: +31-592-31 93 93
info@drweigert.nl

Dr. Weigert UK Ltd
Newmarket Suite 17
Keys Business Village
Keys Park Road
Cannock, WS12 2HA
Tel.: +44-1543478 817
enquiry@drweigert.com

Dr. Weigert Polska Sp. z o.o.
ul. Wybrzeże Gdynskie 6A
01-531 Warszawa
Tel.: +48-22-616 02 23/31
office_pl@drweigert.com

Dr. Weigert España
Avda. de Europa 19, 3A
Parque Empresarial La Moraleja
28108 Madrid
Tel.: +34-91 793 2148
info.spain@drweigert.com

Dr. Weigert France SAS
22 avenue des Nations
Immeuble le Raspail
BP 88035 Villepinte
95932 Roissy CDG Cedex
Tel.: +33-1-48 67 90 33
info@drweigert.fr

Dr. Weigert Hungária Kft.
Hunyadi János út 16
1117 Budapest
Tel.: +36-1-237 06 04
info@drweigert.hu

GB 1023 / 02-2016 © Dr. Weigert, Hamburg

Dr. Weigert management system – Certified according to DIN EN ISO 9001/13485/14001

www.drweigert.com

Reprocessing of Bedpans
and Urine Bottles.

neodisher® SBN plus – Special Rinse Aid for Bedpans.

neodisher® SBN plus

Effective special rinse aid for the reprocessing of bedpans and urine bottles

- With fresh lime blossom fragrance against unpleasant odours
- For washer disinfectors with steam disinfection



5 L canister
4792 35

Fields of application:

- For water softening and rinsing when reprocessing bedpans and urine bottles in bedpan washer disinfectors with steam disinfection
- For optimized, fast and stain-free drying
- Suitable for acid-compatible bedpan washers
- For removing stubborn soiling we recommend the additional use of **neodisher® SBR extra**, a special detergent for bedpans and urine bottles

Disinfection with hot steam:

- With bedpan washers using a thermal reprocessing process the disinfecting step is effected with a hot water-steam mixture
- The steam generator is the heart of the bedpan washer and must be protected against the formation of lime which may cause longer heating times, higher power consumption, overheating and thus could lead to irreparable damage to the heating elements

Our Competence – Your Benefit:

- acidic special rinse aid based on organic acids
- with special softening components for effectively complexing water hardness
- with pleasant fresh lime blossom fragrance
- with surfactants in a low-foam formula



- effectively prevents the formation of lime in the steam generator, in the flushing chamber and on the item to be reprocessed
- suitable for all degrees of water hardness
- reliably covers unpleasant odours
- for optimized, fast and stain-free drying without lime deposits



Pay attention to acid compatibility of dosing systems and steam generators.
For further information on products and applications, please refer to the respective product data sheet.

neodisher® Cleaning – Water Softening & Rinsing.

neodisher® SBR extra

Material-compatible detergent, even removes stubborn residues

- Improves the cleaning performance, removes particularly stubborn residues (such as lime soap stains and residues from cytostatics in faeces)
- Suitable for bedpan washers equipped with a respective dosing pump



200 kg drum
4007 13



5 L canister
4007 35

neodisher® SBK

Combined softening and rinsing agent based on organic acids

- Prevents the formation of lime in the steam generator, pipes and jets
- For fast and optimum stain-free drying
- Suitable for acid-compatible bedpan washers



200 kg drum
4006 13



5 L canister
4006 35

neodisher® SBX

Mildly-alkaline softening agent, complexes water hardness

- Prevents the formation of lime deposits in the steam generator
- Suitable for all bedpan washers



5 L canister
4015 35

