



DR. WEIGERT
Systematic Hygiene

Dental



Our
**Product
Range.**

Cleaning and Disinfection of
Dental Instruments.



neodisher® –
Our Experience for Your
Safety.

For over 30 Years

... neodisher® has been guide to the technical development of instrument hygiene and has set new standards.

neodisher® concepts guarantee first-class results and the longevity of high-quality materials.

Our specific expertise enables us to find solutions to even the most difficult challenges in areas where hygiene plays a crucial role. Put your trust in our expertise.

Expertise from a Single Source.

Reproducible and validated cleaning and disinfection processes are essential for maintaining a safe level of hygiene. Make the most of our expertise: the specialist knowledge of our neodisher® advisers is a strong argument for choosing Dr. Weigert. We are happy to advise you on all issues relating to cleaning and disinfecting medical devices and accessories.

In product and process development neodisher® is always one step ahead. A close co-operation with manufacturers of instruments and other medical devices is the basis for developing innovative process chemicals and designing new processes. Tailor-made practical solutions – for safety and value preservation!

Use Dr. Weigert's know-how to implement hygiene processes in such a way as to put efficiency, safety and cost-effectiveness in optimum relation. Trust in our experience and our service!

Systematic Service for First-Class Results.

Our systematic service provides ongoing support so that you can continually achieve first-class results. We can quickly find answers to all your questions by drawing on our many decades of experience, exceptional expert knowledge and close partnership with the manufacturers of medical devices as well as of washer disinfectors. Your neodisher® specialist advisor is your contact for all questions – providing reliable and professional answers to issues relating to the safe reprocessing of medical devices.



neodisher® Product Range

Manual Reprocessing of Dental Instruments	4
Automated Reprocessing of Dental Instruments:	
Cleaning	5
Neutralisation	6
Rinsing	7
Care of Dental Instruments	7
Surface Cleaning and Surface Disinfection	8
Hand Hygiene	9
Application and Dosing Aids	10
Application Recommendations	11





Manual Reprocessing of Dental Instruments

Correct Reprocessing

... begins with the professional pre-treatment of medical devices – preferably straight after use to prevent the cleaning process from being affected by dried-on residues.

Manual Cleaning and Disinfection	neodisher® Septo PreClean	neodisher® LM 2	neodisher® Septo Plus	neodisher® Septo Fin
cleaning	●	●	●	
disinfecting	●		●	●*
aldehyde-free	●	●	●	
stubborn residues		●		
wet disposal	●			

* also suitable for the final disinfection of semi-critical medical devices

neodisher® Septo PreClean

Detergent and disinfectant – liquid concentrate –

Based on amines – bactericidal, yeasticidal, active against enveloped viruses (1.0 %, 15 min, 20 °C)

For the manual cleaning and disinfection prior to automated reprocessing as well as for wet disposal over longer periods of time, e.g. over night or over the weekend. Very good cleaning performance and excellent material compatibility; no protein fixation. Free of aldehydes and quaternary ammonium compounds. VAH¹ listed and included in the IHO² list of virucidal products.



Canister	5 L	4503 35	Bottles	10 x 1 L	4503 45
----------	-----	---------	---------	----------	---------

neodisher® LM 2

Mildly alkaline detergent – liquid concentrate –

For the cleaning of dental instruments in immersion baths and ultrasonic baths. Supports the removal of stubborn residues such as dental cement. Suitable for stainless steel or chromium-plated instruments. Low-foaming.



Canister	10 L	4403 30	Bottles	10 x 1 L	4403 45
----------	------	---------	---------	----------	---------

neodisher® Septo Plus

Detergent and disinfectant – liquid concentrate –

Based on quaternary ammonium compounds – bactericidal, yeasticidal, active against enveloped viruses (1.0 %, 15 min, 20 °C)

For the manual cleaning and disinfection as well as for the disinfecting pre-cleaning in immersion baths and ultrasonic baths prior to automated reprocessing. Excellent cleaning performance, does not fixate proteins. Excellent material compatibility. Free of aldehydes and amines. Compatible with disinfectants containing aldehyde and peracetic acid. VAH¹ listed and included in the IHO² list of virucidal products.



Canister	5 L	4209 35	Bottles	10 x 1 L	4209 45
----------	-----	---------	---------	----------	---------

neodisher® Septo Fin

Disinfectant – liquid concentrate –

Based on aldehydes – bactericidal, fungicidal, mycobactericidal, virucidal, sporicidal (C. difficile) (5.0%, 60 min, 20 °C)

For use in immersion baths or ultrasonic baths. Free of formaldehyde and glyoxal. Excellent material compatibility. VAH¹ listed, included in the IHO² list of virucidal products, ÖGHMP³ listed.



Canister	5 L	4207 35	Bottles	6 x 2 L	4207 42
----------	-----	---------	---------	---------	---------

¹ Verbund für Angewandte Hygiene [Association of Applied Hygiene]

² Industrieverband Hygiene und Oberflächenschutz [Industry Association Surface Protection and Hygiene, Germany]

³ Österreichische Gesellschaft für Hygiene, Mikrobiologie und Präventivmedizin [Austrian Society of Hygiene, Microbiology and Preventive Medicine]

Automated Reprocessing of Dental Instruments



Automated Reprocessing

... is done using alkaline detergents, which – depending on product type – can also contain surfactants and enzymes. Surfactants have emulsifying and dispersing properties, enabling substances that are insoluble in water to be emulsified in water and suspended in the solution. Enzymes enable the breakdown of organic residues.

Automated Cleaning	neodisher® MediClean Dental	neodisher® FA Dental	neodisher® MA Dental
alkaline	●	●	●
with surfactants	●		●
enzymatic	●		
liquid	●	●	
excellent material compatibility	●		
for Miele and Melag washer disinfectors	●	●	●

neodisher® MediClean Dental

Mildly alkaline, enzymatic detergent – liquid concentrate –

Effectively removes residues of blood and protein. Very gentle on material. Also suitable for the manual reprocessing of transmission instruments. Recommended for Miele and Melag washer disinfectors.

Canister 5 L 4053 33



neodisher® FA Dental

Alkaline detergent – liquid concentrate – Free of surfactants and oxidative additives

Good cleaning performance with excellent material compatibility. Recommended for Miele and Melag washer disinfectors.

Canister 5 L 4124 33



neodisher® MA Dental

Alkaline detergent – powder – Free of oxidative additives

Good cleaning performance with excellent material compatibility. Recommended for Miele and Melag washer disinfectors.

Bucket 10 kg 4113 76 Bottles 12 x 1 kg 4113 87





Automated Reprocessing of Dental Instruments

Our Recommendations

... for the neutralisation of alkaline residues in the automated cleaning of dental instruments in washer disinfectors

Neutralisation	neodisher® N Dental	neodisher® Z Dental
acidic	●	●
free of surfactants	●	●
free of phosphates		●
for sensitive materials		●

neodisher® N Dental

Neutralising agent – liquid concentrate – Based on phosphoric acid. Free of surfactants

Neutralises carried over alkali residues and removes discolourations, lime residues and rust during automated reprocessing of dental instruments. Also suitable for thorough cleaning of dental instruments in dip baths. Recommended for Miele and Melag washer disinfectors.



Canister 5 L 4118 33

Bottles 12 x 1 L 4118 87

neodisher® Z Dental

Neutralising agent – liquid concentrate – Based on organic acids. Free of surfactants and phosphates

Neutralises carried over alkali residues during automated reprocessing of dental instruments. Also suitable for the automated reprocessing of transmission instruments. Recommended for Miele and Melag washer disinfectors.



Canister 5 L 4055 33

Bottles 12 x 1 L 4055 87



Automated Reprocessing & Care of Dental Instruments



Rinsing and Instrument Care	neodisher® MediKlar Dental	neodisher® Mielclear	neodisher® IP Spray
Rinse aid	●	●	
Instrument lubricant			●
for automated use	●	●	
for manual use			●

Rinse aids

... enable an improved, stainfree drying of dental instruments and utensils when reprocessing automatically.

neodisher® MediKlar Dental

Rinse aid – Liquid concentrate –

With special wetting agents for an improved drying of dental instruments incl. transmission instruments and utensils when reprocessing automatically. Also suitable for difficult-to-wet materials e.g. synthetic materials. Recommended for Miele and Melag washer disinfectors.



Bottles 12 x 1 L 4776 49

neodisher® Mielclear

Rinse aid – Liquid concentrate –

For the rapid and stainfree drying of dental instruments and utensils when reprocessing automatically. With corrosion inhibitor. Recommended for Miele washer disinfectors.



Bottles 12 x 1 L 4119 49

neodisher® IP Spray

Lubricant for instruments – spray – Based on highly purified white oil

For the manual lubrication of hinged instruments after automated reprocessing. Does not lead to encrustations, does not influence a following steam or hot air sterilisation. Toxicologically harmless.



Spray can 12 x 400 ml 4304 90

Surface Cleaning and Surface Disinfection

neoform MED AF*

Aldehyde-free disinfecting cleaner - liquid concentrate - bactericidal, yeasticidal, active against enveloped viruses

For cleaning and disinfecting medical devices and medical inventory. Application: wiping method. VAH¹ listed and included in the IHO² list of virucidal products.

Canister 5 L 4056 35



neoform MED FF*

Disinfecting cleaner – liquid concentrate – bactericidal, tuberculocidal, yeasticidal, virucidal, active against C. difficile

For cleaning and disinfecting medical devices and medical inventory. Application: wiping method. VAH¹ listed and included in the IHO² list of virucidal products.

Canister 5 L 4057 35



neoform MED rapid

Rapid disinfectant – ready-to-use solution – bactericidal, yeasticidal, tuberculocidal, active against enveloped viruses, active against Norovirus

Alcohol-based solution for disinfecting medical devices. Application: wiping method and spraying devices. VAH¹ listed and included in the IHO² list of virucidal products.

Bottles 12 x 0,75 L 4070 47



Hand Hygiene



triformin® wash

Washing lotion

For hand and body washing. Excellent skin tolerance, confirmed by dermatological certificates. Especially suitable for normal, dry and sensitive skin. Based on natural plant-based raw materials, with aloe vera and natural bee honey. Free of soap and alkali. Dosing via triformin® dispenser 1.0 and triformin® dispenser 0.5.



Canister	5 L	3173 35	Bottles	20 x 500 ml	3173 63
Bottles	10 x 1 L	3173 48			

triformin® care

Skin care cream

For skin care after work and during breaks. Excellent skin tolerance, confirmed in dermatological certificates. Free of perfumes. Dosing of 500 ml bottle via triformin® dispenser 0.5.



Bottles	6 x 500 ml	3105 39	Tubes	12 x 100 ml	3105 15
---------	------------	---------	-------	-------------	---------

triformin® protect

Skin protection cream

Protects the skin from drying out. Long-lasting protection and intensive care. Excellent skin tolerance, confirmed by dermatological certificates. Free of perfumes. Dosing of 500 ml bottles via triformin® dispenser 0.5.



Bottles	6 x 500 ml	3102 39	Tubes	12 x 100 ml	3102 15
---------	------------	---------	-------	-------------	---------

Application and Dosing Aids

Dosing and Application Ancillaries for Canisters and Drums

Plastic key for the easy opening of canisters and drums
 Canister pump for 5 L and 10 L canisters (30 ml stroke)
 Outlet tap for 5 L and 10 L canisters
 Spraying attachment for neoform MED rapid

Art.-No. 9810 07
 Art.-No. 9807 01
 Art.-No. 9800 04
 Art.-No. 9500 07



Dosing and Measuring Beakers

500 ml measuring beaker, graduated
 1,000 ml measuring beaker, graduated

Art.-No. 9801 06
 Art.-No. 9801 00



Personal Protective Clothing

Protective gloves
 Protective goggles
 Eye-wash bottle

Art.-No. 9880 05
 Art.-No. 9880 06
 Art.-No. 9880 07



Disinfection Tubs

3 l disinfection tub (with sieve insert and lid)
 10 l disinfection tub (with sieve insert and lid)
 30 l disinfection tub with drain tap (with sieve insert and lid)

Art.-No. 9611 00
 Art.-No. 9611 20
 Art.-No. 9611 40



Dosing Dispensers and Hand Pumps for triformin® Products

Dispenser
 triformin® dispenser 1.0 gold* for dosing from 1 L bottles
 triformin® dispenser 0.5 gold* for dosing from 500 ml bottles
 triformin® dispenser 1.0 silver** for dosing from 1 L bottles
 triformin® dispenser 0.5 silver** for dosing from 500 ml bottles
 triformin® dispenser 1.0 basic*** for dosing from 1 L bottles
 triformin® dispenser 0.5 basic*** for dosing from 500 ml bottles

Art.-No. 9810 51
 Art.-No. 9810 21
 Art.-No. 9810 52
 Art.-No. 9810 22
 Art.-No. 9810 53
 Art.-No. 9810 23

Hand pumps
 for 1 l bottles triformin® wash
 for 0.5 l bottles triformin® wash
 for 0.5 l bottles triformin® care and triformin® protect

Art.-No. 9802 00
 Art.-No. 9802 03
 Art.-No. 9802 02

Wall bracket
 for 1 l bottles triformin® wash
 for 0.5 l bottles triformin® wash, triformin® care, triformin® protect

Art.-No. 9810 56
 Art.-No. 9810 59



* washer-proof and autoclavable; stainless steel pump easily removable from the front
 ** washer-proof and autoclavable, with stainless steel pump
 *** with plastic pump





Application Recommendations

Automated Reprocessing	Dental instruments and utensils	Transmission instruments
Cleaning		
neodisher® MediClean Dental	●	●
neodisher® FA Dental	●	
neodisher® MA Dental	●	
Neutralisation		
neodisher® Z Dental	●	●
neodisher® N Dental	●	
Rinsing		
neodisher® MediKlar Dental	●	●
neodisher® Mielclear	●	●

Our Recommendations

The lefthand table lists optimum combinations of process chemicals depending on the instruments to be reprocessed. We recommend a combination of neodisher® MediClean Dental, neodisher® Z Dental and neodisher® MediKlar Dental, which is equally suitable for reprocessing dental instruments, dentist's utensils (trays, rinsing beakers, kidney dishes) as well as transmission instruments thanks to its excellent material compatibility.

Hand Hygiene	triformin® wash	triformin® care	triformin® protect
cleaning	●		
care		●	
protection			●

Manual Cleaning, Disinfection, Care	neodisher® Septo PreClean	neodisher® LM 2	neodisher® Septo Plus	neodisher® Septo Fin	neodisher® N Dental	neodisher® IP Spray	neoform MED AF**	neoform MED FF**	neoform MED rapid
pre-treatment of dental instruments	●		●						
manual cleaning of dental instruments	●	●	●						
manual disinfection of dental instruments	●		●	● *					
thorough cleaning of dental instruments					●				
care of dental instruments						●			
surface cleaning							●	●	
surface disinfection							●	●	●

* also suitable for the final disinfection of semi-critical medical devices
 ** Use disinfectants safely. Always read the label and product information before use.



DR. WEIGERT
Systematic Hygiene

Always on Hand
to Help with Hygiene –
**in Europe and Around
the World.**



Chemische Fabrik
Dr. Weigert GmbH & Co. KG
Mühlenhagen 85
20539 Hamburg
Tel.: +49-40-7 89 60-0
info@drweigert.de

Dr. Weigert
Handelsgesellschaft m.b.H.
Wienerbergstraße 11/12 a
1100 Wien
Tel.: +43-1-99460 6275
info@drweigert.at

Dr. Weigert (Schweiz) AG
Baarerstrasse 14
6300 Zug
Tel.: +41-41 561 32 90
info@drweigert.ch

Dr. Weigert Nederland BV
Balkendwarsweg 3
9405 PT Assen
Tel.: +31-592-31 93 93
info@drweigert.nl

Dr. Weigert UK Ltd
Newmarket Suite 17
Keys Business Village
Keys Park Road
Cannock, WS12 2HA
Tel.: +44-1543 471 222
enquiry@drweigert.com

Dr. Weigert Polska Sp. z o.o.
ul. Wybrzeże Gdynskie 6D
01-531 Warszawa
Tel.: +48-22-616 02 23/31
office_pl@drweigert.com

Dr. Weigert España
Avda. de Europa 19, 3A
Parque Empresarial La Moraleja
28108 Madrid
Tel.: +34-91 793 2148
info.spain@drweigert.com

Dr. Weigert France SAS
22 avenue des Nations
Immeuble le Raspail
BP 88035 Villepinte
95932 Roissy CDG Cedex
Tel.: +33-1-48 67 90 33
info@drweigert.fr

Dr. Weigert Hungária Kft.
Hunyadi János út 16
1117 Budapest
Tel.: +36-1-237 06 04
Info@drweigert.hu

Dr. Weigert management system – Certified according to DIN EN ISO 9001/13485/14001

www.drweigert.com